

Grace[®] Silica Drug Delivery

Increasing Drug Solubility with Solid Dispersions of SilSol[™] 6035 Silica Formulated Into Sublingual Tablets

Introduction

Improving the solubility of BCS class 2 drugs is a continuing challenge for the pharmaceutical formulator. Improving solubility could allow for formulations using significantly less active and open the door to new dosage forms and avenues for lifecycle extensions. Creating solid dispersions are a promising way to improve drug solubility. Grace SilSol[™]6035 has a specifically engineered surface and pore distribution to enable stable amorphous dispersions by solvent deposition. Poorly soluble Nifedipine was loaded on SilSol[™] 6035 silica to make a solid dispersion with enhanced water solubility. Nifedipine is stabilized in its amorphous form by confinement effects of the engineered pore size of the silica. This solid dispersion greatly improves the solubility of Nifedipine and, since it is a free flowing powder, is easily formulated into solid oral dosage forms such as tablets.

Materials and Methods

1. Nifedipine was loaded on SilSol[™] 6035 using solvent method of drug loading: saturated solution of Nifedipine was prepared in acetone and adsorbed on SilSol[™] 6035 drop wise with continuous stirring and drug was allowed to penetrate in Silica pores. Acetone was evaporated. Dissolution study was conducted to confirm the successful drug loading and increased water solubility.

Release of Nifedipine from SilSol[™] 6035 dispersion

The release study was conducted in water as a dissolution media using USP apparatus II. SilSol[™] helps to get 100% of Nifedipine released in first 5min while less than 10% was released from the pure drug powder.




Figure 1: Dissolution of Nifedipine from SilSol[™] 6035 and alone drug

Grace[®] Silica Drug Delivery



- The tablets of obtained solid dispersion was compressed using formula mentioned below:

| Ingredients | mg/tab |
|--|--------|
| Nifedipine Dispersion (10% drug on SilSol [™] 6035) | 100 |
| Mannitol | 140.85 |
| Cross povidone | 37.8 |
| Poly vinyl acetate | 7.85 |
| Sodium saccharine | 3 |
| Syloid 244 FP | 3 |
| Magnesium stearate | 1.5 |
| Citric acid | 6 |
| Total | 300 |

Table 1: Formulation of Nifedipine sublingual tablets using SilSol[™] 6035 and Syloid[®] 244FP

- Syloid[®] 244FP was used as glidant which also helps to increase saliva diffusion in tablets.
- Lubricated powder was subjected to compression by using Eliza Press 200 tablet press.
- In process quality control parameters were evaluated for the manufactured tablets.

Results

The prepared tablets were evaluated for various parameters and observations of the same are reported in table below:

| Parameter | Mean ±SD |
|---------------------------------------|----------|
| Weight variation (mg) | 299.35 |
| Hardness (Kg/cm ²) | 6.52 |
| Thickness (mm) | 5.67 |
| Friability (%) | 0.07 |
| Wetting time (sec) | 22.21 |
| <i>In vitro</i> dispersion time (sec) | 41.32 |
| Water absorption ratio (%) | 53.22 |

Table 2: Observations of Nifedipine sublingual tablets

Conclusions

SilSol[™] 6035 helps to increase the solubility of poorly soluble drugs like Nifedipine. Drug loaded SilSol[™] can easily compressed into sublingual tablets using Syloid[®] 244FP as glidant. The hydrophilic surface and Porous nature of SilSol[™] 6035 help to increase drug release for immediate absorption. The platform technology can be effective to increase bioavailability of poorly soluble drugs from BCS II and IV class by creating solid dispersions and formulating into solid dosage forms.

GRACE[®] is a trademarks, registered in the United States and/or other countries of W. R. Grace & Co.-Conn. SILSOL[™] is a trademark of W. R. Grace & Co.-Conn. This trademark list has been compiled using available published information as of the publication date of this brochure and may not accurately represent current ownership/or status. Grace Materials Technologies is a business segment of W. R. Grace & Co.-Conn., which includes products formerly sold under the Alltech and Grace Davison brands. The information presented herein is derived from our testing and experience. It is offered for your consideration and verification. Since operating conditions vary significantly, and are not under our control, we disclaim all warranties on the results that may be obtained from the use of our products. W. R. Grace & Co.-Conn. and its subsidiaries can not be held responsible for any damage or injury occurring as a result of improper installation or use of its products. Grace reserves the right to change prices and/or specifications without prior notification. © Copyright 2015 W. R. Grace & Co. - Conn. All rights reserved.

AP024