



Moisture Activated Dry Granulation (MADG) with Syloid® XDP Silica

Introduction

Moisture activated dry granulation (MADG) has been presented to pharmaceutical formulators as an innovative way to modify wet granulation. It provides a way to address common problems such as endpoint detection, drying and milling ¹. This technique offers an efficient, cost-effective granulating process utilizing a limited amount of water and does not require a drying step ². Because this process depends on the addition of a moisture absorber to facilitate even moisture distribution and drying of the granulated powder, careful selection of excipients is critical for the effective utilization of this technique ³. Syloid® XDP mesoporous silica excipient is an effective moisture absorber for this technique. The high surface area and pore volume coupled with an active silica surface makes Syloid® XDP silica a highly efficient adsorbent for moisture, effectively “locking” it away in its pores to create a stable, dry powder. In addition, the large particle size of Syloid® XDP silica helps promote a blend with improved flow characteristics for easier processing.

MADG Formulation

A MADG formulation was developed by using Paracetamol as the API of interest with the final dose content of 200mg per tablet. Initially API, fillers and dry binder were added into a rapid mixer-granulator followed by the addition of 5% water spray onto the dry mixture, allowing the process to form granules by activation of dry binder upon contact with the water spray. After 3-5 minutes of granulation, Syloid® XDP silica is added to distribute the moisture evenly and dry the granulated powder. Once the granulated powder is dry then the remaining excipients, such as disintegrant, glidant and lubricant are added. The final lubricated blend was subjected for compression study to evaluate the in-process quality control parameters. For the comparison purpose the same composition was used for direct compression without any addition of water.

Formula


| Ingredients | Description | Composition % w/w |
|----------------------------------|--------------|-------------------|
| Paracetamol | API | 40 |
| Di-Calcium Phosphate (anhydrous) | Filler | 12.5 |
| Poly Vinyl Pyrrolidone K-30 | Dry Binder | 5 |
| Micro crystalline cellulose 102 | Filler | 33.5 |
| Syloid® XDP3050 silica | Adsorbent | 2.5 |
| Ac-Di-Sol | Disintegrant | 5 |
| Syloid® 244FP silica | Glidant | 1 |
| Magnesium Stearate | Lubricant | 0.5 |

Flow Parameters

| Flow Parameters | MADG | Direct Compression |
|---------------------------|-------|--------------------|
| Angle of Repose° | 35.7 | 44.35 |
| Bulk Density (g/mL) | 0.417 | 0.448 |
| Tapped Density (g/mL) | 0.517 | 0.625 |
| Compressibility Index (%) | 19.44 | 28.35 |
| Hausners Ratio | 1.241 | 1.396 |

Tabletting Parameters

| Tabletting Parameters | MADG | Direct Compression |
|-----------------------|-------|--------------------|
| Average Weight (mg) | 501.3 | 501.7 |
| Hardness (N) | 125.2 | 50.2 |
| Thickness (mm) | 4.68 | 5.13 |
| Disintegration (min) | 0.50 | 0.50 |
| Friability (%) | 0.27 | 0.35 |



Conclusion

- MADG process is advantageous over conventional wet granulation in terms of reduced process time.
- Syloid® XDP silica acts as an effective adsorbent which eliminates the drying step.
- Syloid® XDP silica helps to distribute drug uniformly in the formulation and generates uniform sized granules.
- Tablet blend can be generated in single piece equipment (rapid mixer granulator).
- Flow characteristics of the powder blend is improved vs direct blending.
- Tablet hardness is improved vs direct compression.
- Useful for moisture sensitive drugs where wet granulation is not possible.

References

- 1 Ullah et al., Pharmaceutical Technology, Vol.II, 48-54, 1987
- 2 Ullah et al., Pharmaceutical Technology, PharmTec, March 7, 2011
- 3 Ullah et al., Pharmaceutical Technology, November 2009



Global Headquarters
7500 Grace Drive
Columbia, MD 21044 USA

Americas
Deerfield, Illinois
Tel: +1 847.948.8600

Europe
Lokeren, Belgium
Tel: +32 (0)9.340.65.65

India
Pune
Tel: +91 20.6644.9900

China
Shanghai
Tel: +86 21.3158.2888

Australia
Melbourne
Tel: +61 3.8401.6369

The information presented herein is derived from our testing and experience. It is offered for your consideration and verification. Since operating conditions vary significantly, and are not under our control, we disclaim all warranties on the results that may be obtained from the use of our products. W. R. Grace & Co.-Conn. and its subsidiaries can not be held responsible for any damage or injury occurring as a result of improper installation or use of its products. Grace reserves the right to change prices and/or specifications without prior notification. GRACE® and SYLOID® are trademarks, registered in the United States and/or other countries, of W. R. Grace & Co.-Conn. This trademark list has been compiled using available published information as of the publication date of this brochure and may not accurately reflect current trademark ownership or status. Grace Discovery Sciences is a product group of W. R. Grace & Co.-Conn. Alltech Associates, Inc. is a wholly owned subsidiary of W. R. Grace & Co.-Conn. Grace Materials Technologies is a business segment of W. R. Grace & Co.-Conn. which includes products formerly sold under the Alltech and Grace Davison brands. © Copyright 2015 W. R. Grace & Co.-Conn. All rights reserved.